

PalmTop™ docking station specifications

Ventilator docking interface

Ventilator communication port	The PalmTop™ ventilator docking station will establish communication with a ventilator within five seconds, when docked	
	VOXP protocol	
Nurse call port	Mechanical	4-position RJ type female connector
	Electrical	Closed contact resistance ≤ 1 Ohm, electrically isolated Open contact resistance ≥ 1 M.Ohm, electrically isolated

Power

External AC power connection	Input voltage range	100 to 240 VAC, 50 to 60 Hz
	Input power connector	IEC 320 (with connector clamp)
PalmTop Monitor DC power connection	Output voltage	14 to 16 VDC
	Output current	up to 1.4 A





Physical dimensions (width x height x depth)


Size	5 5/8" w x 7 1/4" h x 12" d (14.2 cm x 18.4 cm x 30.4 cm)
Weight	5 lbs (2.3 kg)

Environmental

Temperature	Storage	-20° C to 60° C (-4° F to 140° F)
	Operating	5° C to 40° C (41° F to 104° F)
Humidity	Storage	Up to 95% relative, non-condensing
	Operating	15% to 95% relative, non-condensing
Altitude	Storage	50,000' maximum above sea level
	Operating	68 to 106 kPa

 **WARNING:** U.S. Federal Law restricts this device to sale by or on the order of a physician.

 CareFusion
22745 Savi Ranch Parkway
Yorba Linda, CA 92887
800.231.2466 toll-free
714.283.2228 tel
714.283.8493 fax

 CareFusion Germany 234 GmbH
Leibnizstrasse 7
97204 Hoechberg
Germany
+49 931 4972-0 tel
+49 931 4972-423 fax

carefusion.com

© 2011 CareFusion Corporation or one of its subsidiaries. All rights reserved. PalmTop is a trademark or registered trademark of CareFusion Corporation or one of its subsidiaries. RC2888 (0711/4000) L3427 Rev. A

